

TheraFace Hot & Cold Rings Product Description / PDP Copy

Hot & Cold Rings:

Powered by Cryothermal Technology™

The TheraFace Hot and Cold Rings, which attach to your TheraFace PRO device, feature a combined 6 varying temperature settings to better help your skin absorb creams, help reduce facial puffiness, improve radiance, and much more.

Hot Ring:

- Relaxes and reduces tension
- Decreases muscle spasms
- Increases blood flow

Cold Ring:

- Decreases pain
- Decreases inflammation
- Decreases muscle spasms

PRO tip: The Cold Ring attachment can be used with your favourite percussive attachment for the ultimate treatment to depuff and relax your facial muscles.

Features:

- Compatible with the TheraFace PRO
- Hot Ring Temperature Settings:
 - 35C/95F
 - 39C/102F
 - 43C/109F
- Cold Ring Temperature Settings:
 - 26C/78F
 - 22C/71F
 - 18C/65F
- Available in two colour-ways, White or Black
- Easy to clean for a hygienic experience

What's Included:

- 1 TheraFace Hot Ring
- 1 TheraFace Cold Ring
- Soft Carrying Case