



FitZone

s o l u t i o n s

R a i s i n g t h e b a r

Prepared by FitZone Solutions
COPYRIGHT 2024



Scan & Save – Your Digital Blood Pressure Diary

for all Beurer blood pressure monitors – even without Bluetooth®

Writing down your blood pressure values is a thing of the past! Thanks to the "Scan & Save" function, you can now save measurements from all Beurer blood pressure monitors digitally in the app – even without Bluetooth*. All you need to do is scan the measurements with a smartphone camera. The result is automatically detected and saved in the app "beurer HealthManager Pro".

Just 4 easy steps....

beurer

Your Digital Blood Pressure Diary
in a few simple steps:

1. Download the App
"beurer HealthManager Pro"




Scan & Save
Digital Blood Pressure Diary

beurer

Quick

2. Open the "Scan & Save"
function in the app and scan the measurements




Scan & Save
Digital Blood Pressure Diary

beurer

Easy

3. The measurements are automatically saved
in the app



Scan & Save
Digital Blood Pressure Diary

beurer

Even without *Bluetooth*®.

4. Keep track of your daily measurements
in the app at any time.



Scan & Save
Digital Blood Pressure Diary